



BM-TEST

USER GUIDE

H10798 - Version 1 / 03-2018

1 - Product presentation

The magnetic loop earpiece is a compact device that allows people without hearing aids to listen to a hearing loop system as if they were wearing a hearing aid in the "T" position.

It can be used for:


- Check the sound source of a magnetic loop
- Equipping a person who is hard of hearing and who does not have a hearing aid
- Equipping a person with hearing loss whose prosthesis does not have a "T" position.

2 - Components

Headphone

Volume control knob

LED



Hearing loop listener

Note: The BM-TEST is not a magnetic field meter and does not allow EN60118-4 commissioning.

3 - Operation

The hearing loop listener works with two AAA alkaline batteries.

1. Place the headset on your ears and adjust it for comfortable listening.
2. Hold the hearing loop listener vertically and at head height to mimic the position of a hearing aid.



3. Turn the unit on by turning the volume control slightly clockwise until a comfortable listening sound is obtained. Turn the knob in the opposite direction to decrease the volume if necessary..
4. If the LED does not light, check that the batteries work and are correctly placed as shown inside the battery.
5. If the LED lights up but you can't hear a sound, ask a colleague or user to speak into the microphone position on the loop system or test the input source if the sound is connected. Re-check the hearing loop listener reception.
6. When you are finished using the headset, turn the volume knob counterclockwise until the LED turns off.

Important

Always turn off the unit when it is not in use to extend its life.

Always hold the hearing loop listener at ear height and upright when listening to a loop system to ensure the best possible performance.

4 - Maintenance

To extend the life of the product, perform the following operations regularly:

Clean the unit with a soft, dry cloth or suitable cleaning wipes. Do not use solvent-based cleaners and avoid getting all kinds of liquids into the hearing loop listener.

Turn off the headset when not in use.

5 - Troubleshooting

Symptom	Possible fault	Action
Ignition LED does not light up.	<ol style="list-style-type: none">1. The batteries are empty2. There is an internal fault	<ol style="list-style-type: none">1. Replace the batteries with new ones (AAA alkaline)2. Contact your dealer
Whistle heard through the loop	Interferences	At the top of the battery compartment of the device, you can find a noise filter switch accessible by removing the batteries. Move the switch to the right to turn it on.
No loop sound is heard in the hearing loop listener.	<ol style="list-style-type: none">1. Sound is not loud enough.2. The earphone is not at an appropriate distance to listen to the loop.3. The magnetic loop does not work.4. The headphone has a problem.5. Internal fault.	<ol style="list-style-type: none">1. Turn up the volume.2. Bring the hearing loop listener closer to the loop.3. Use another hearing loop listener or listening system (if possible) to make sure the loop works. If no sound is coming out of the loop, contact your dealer.4. Replace the headphone.5. Follow step 3. If the sound is heard with another device, contact your dealer to change the hearing loop listener.

Because AUDIOPHONY® takes the utmost care in its products to make sure you only get the best possible quality, our products are subjects to modifications without prior notice. That is why technical specifications and the products physical configuration might differ from the illustrations.

Make sure you get the latest news and updates about the AUDIOPHONY® products on www.audiophony.com
AUDIOPHONY® is a trademark of **HITMUSIC SAS** - Parc d'activités Cahors sud - 46230 FONTANES - FRANCE